







Accessories for 2.853.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.853.0010

Qty.	Order no.	Description	
1 pcs	1.853.0010	853 «MCS» – CO₂ Suppressor «MCS» – CO ₂ Suppressor for sequential suppression. Retrofits all existing Metrohm Advanced IC Systems	
1 pcs	6.1801.140	PVC Tubing 3.2 / 6.4 mm / 11 cm Connecting CO ₂ adsorber cartridge and H ₂ O adsorber cartridge of the CO ₂ Suppressor.	
1 pcs	6.1808.190	Adapter tubing nozzle/Luer For connecting the CO ₂ Adsorber Cartridge and H ₂ O Adsorber Cartridge	
1 pcs	6.2027.070	Column holder 25 mm For mounting columns in column chamber of IC instruments.	
1 pcs	6.2128.180	Cable 762/830 - 752/753/754/833 Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit	

1 pcs	6.2143.230	Cable Compact IC - CO₂ Suppressor Connection cable from Compact IC (761 and 861) to CO ₂ Suppressor	
1 pcs	6.2152.020	Power supply / 100-240 V/15 V DC Power supply to 853 MCS Suppressor	
2 pcs	6.2837.010	H₂O adsorption cartridge To CO ₂ Suppressor. Water adsorption cartridge for the aspirated air.	
1 pcs	6.2837.100	CO₂ adsorption cartridge CW Adsorption cartridge for removing the CO ₂ from the air aspirated by the MCS. The cartridge is used without upstream H ₂ O adsorption cartridge.	