



## Accessories for 6.02041.020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.01404.020</b>	<b>Beaker adapter for 200 mL sample beaker</b>	
	Beaker adapter for centering the 200 mL sample beaker in the slide of the Pick&Place modules.	
<b>6.1459.310</b>	<b>Sample beaker, disposable, plastic (PP), 200 mL, 1000 pieces</b>	
	Sample beaker, disposable, plastic (polypropylene, PP), 200 mL, suitable for 6.02041.020 OMNIS Software sample rack, 6.2041.370 und 6.2041.830 sample racks, 1000 pieces.	