





Accessories for 6.02601.040

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces	
Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.		
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces	
Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.		
6.01400.303	Sample beaker (20x) amber glass 120 mL (P&P)	
20 sample beakers 120 mL made of amber glass for the OMNIS Sample Robot Pick&Place		
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces	
Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces		

6.01402.003 Sample beakers, amber glass, 75 mL, 30 pieces

Sample beakers, amber glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces



6.05700.000 Consumable Kit OMNIS Gripper

Four slide sleeves and rubber gripper fingers each for replacement on the gripper finger of the OMNIS Sample Robot.

