






Accessories for 6.05300.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05300.000

Qty.	Order no.	Description	
1 pcs	6.1015.000	Metrosep C Trap 1 - 100/4.0	
		High-capacity cation column, which is used to purify the eluent flow. Even reagents of the highest quality level, e.g. "Ultrapur" or "puriss", can still contain minimum amounts of cation contamination. These are reliably held back by the Metrosep C Trap 1 - 100/4.0.	
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents	
		For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
1 pcs	6.1608.070	Eluent bottle / 2 L / GL 45	
1 pcs	6.1618.020	Thread adapter S 40 to GL 45	
		Adapter for bottles with S40 thread (Merck...) to GL45.	
1 pcs	6.1624.000	Adapter SGJ 14 for 6.1619.XXX Adsorber tube	
		Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.	

1 pcs **6.1805.120** **FEP tubing / M6 / 100 cm**
With light and kink protection



1 pcs **6.1829.020** **FEP aspiration tube M6, 0.5 m**
For Dosing Unit



1 pcs **6.1831.030** **PEEK capillary 0.75 mm i.d. / 3 m**
For all IC components.



1 pcs **6.1831.120** **PEEK capillary 0.25 mm i.d., 45 cm**



1 pcs **6.1834.090** **Connecting tubing for degasser pump, 125 mm**
For Professional IC Instruments



1 pcs **6.2053.030** **Mounting ribbon to 6.2744.340, 25 x 50 mm**
Mounting ribbon to Equipment set Dosino gradient



1 pcs **6.2057.210** **Holder for Dosino to IC instruments**
Mounting part for Dosinos at the IC instrument housing.



1 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs **6.2744.080** **M6 thread / UNF 10/32 adapter**

With M6 outer thread. For IC instruments.



1 pcs **6.2744.340** **6-port adapter 4 × UNF 10/32, 2 × UNF 1/4-28**

For connecting up to five Dosinos with various solutions to one IC pump. For Dosino gradient applications.



1 pcs 6.2744.350 PEEK pressure screw 14 mm, 3 pcs.



With UNF 10/32 connector. For the connection of PEEK capillaries in the Absorber Module.
