

Accessories for 6.1432.210

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1458.020	Titration head / 4 x M10 connectors	
	With 4 M10 screw nipples, including stoppers.	
6.1458.030	KF titration head / 4 x M10 connectors	
	With 4 M10 screw nipples, including stoppers.	
6.1458.070	Titration head / 2 x SGJ 14	
	Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.	
6.1909.030	Stirring propeller / 104 mm	
	Propeller stirrer with PP shaft, PVDF tip. Length from the SGJ: 104 mm. For KF applications with 722, 802 Rod Stirrer	
6.2037.040	Foil holder	
	Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer	

6.2041.340 **Sample rack 24 x 75 mL**

Sample rack for 24 x 75 mL samples, one-row.



6.2041.350 **Sample rack 48 x 75 mL**

Sample rack for 48 x 75 mL samples, two-row



6.2041.800 **Sample rack 100 x 75 mL**

Sample rack for 100 x 75 mL samples, multi-row



6.2041.900 **Sample rack insert 54 x 75 mL for water bath**

Insert for waterbath rack 54 x 75 mL



6.2820.000 **Al foil**

1000 circular sheets of 80 mm diameter. Thickness: 0.010 mm. Used together with 6.2037.030 Foil holder for 6.1432.210 Sample beaker. For 674 KF Sample Changer

