

## Accessories for 6.5339.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5339.010

Qty.	Order no.	Description	
1 pcs	6.0343.100	<b>Separate Pt Rod Electrode</b>  Platinum auxiliary electrode for voltammetry and CVS.	
1 pcs	6.0728.130	<b>LL-Ag/AgCl reference electrode for CVS</b>  Inner reference system with electrolyte c(KCl) = 3 mol/L. For use with electrolyte vessel made of plastic 6.1245.010.	
1 pcs	6.1204.510	<b>Driving axle for rotating disk electrode (RDE)</b>  Driving axle for rotating disk electrode (RDE) with titanium axle for Professional VA/CSV instruments, M3 thread.	
1 pcs	6.1204.610	<b>2 mm platinum electrode tip for CVS</b>  Metal electrode tip made of platinum (Pt), diameter of electrode disk 2 mm ±0.05 mm, polished surface, shaft made of glass, M3 thread connection. Concentricity error: < 0.3 mm. Temperature range: 0 - 50 °C.	
1 pcs	6.1245.010	<b>Electrolyte vessel made of plastic</b>  Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.	

---

1 pcs    6.1415.210    Titration vessel / 10-90 mL



---

1 pcs    6.1415.250    Titration vessel / 50-150 mL



---

1 pcs    6.2308.020    **Electrolyte 3 mol/L KCl (250 mL)**  
Electrolyte solution  $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)



---

1 pcs    6.2310.010    **Electrolyte KNO<sub>3</sub>-1mol/L 250 mL**  
Electrolyte solution KNO<sub>3</sub> 1 M (reference electrolyte for combined silver electrodes)

