

## Accessories for 6.5330.150

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

### Scope of delivery 6.5330.150

Qty.	Order no.	Description	
1 pcs	6.1014.000	<b>Metrosep A Trap 1 - 100/4.0</b>	 <p>High-capacity anion column, which is used to purify the eluent flow. Even chemicals of the highest quality level, e.g. "Ultrapur" or "puriss", can still contain minimal anion contaminations. These are reliably held back by the Metrosep A Trap 1 - 100/4.0. Should these anionic contaminants elute at a later time, there will be practically no further noticeable influence on the baseline, due to the enormous peak widening.</p> <p>This column is primarily used with gradient applications.</p>
1 pcs	6.1602.160	<b>Bottle attachment GL 45 for eluents and reagents</b>	 <p>For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.</p>
1 pcs	6.1608.070	<b>Eluent bottle / 2 L / GL 45</b>	
1 pcs	6.1618.020	<b>Thread adapter S 40 to GL 45</b>	 <p>Adapter for bottles with S40 thread (Merck...) to GL45.</p>

---

**1 pcs**     **6.1624.000**     **Adapter SGJ 14 for 6.1619.XXX Adsorber tube**

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.



---

**1 pcs**     **6.1805.120**     **FEP tubing / M6 / 100 cm**

With light and kink protection



---

**1 pcs**     **6.1829.020**     **FEP aspiration tube M6, 0.5 m**

For Dosing Unit



---

**1 pcs**     **6.1831.030**     **PEEK capillary 0.75 mm i.d. / 3 m**

For all IC components.



---

**1 pcs**      **6.1831.120**      **PEEK capillary 0.25 mm i.d., 45 cm**



---

**1 pcs**      **6.1834.090**      **Connecting tubing for degasser pump, 125 mm**  
For Professional IC Instruments



---

**1 pcs**      **6.2053.030**      **Mounting ribbon to 6.2744.340, 25 x 50 mm**  
Mounting ribbon to Equipment set Dosino gradient



---

**1 pcs**      **6.2057.210**      **Holder for Dosino to IC instruments**  
Mounting part for Dosinos at the IC instrument housing.



---

**1 pcs**      **6.2744.070**      **Pressure screw short**  
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



---

**1 pcs**      **6.2744.080**      **M6 thread / UNF 10/32 adapter**

With M6 outer thread. For IC instruments.



---

**1 pcs**      **6.2744.340**      **6-port adapter 4 × UNF 10/32, 2 × UNF 1/4-28**

For connecting up to five Dosinos with various solutions to one IC pump. For Dosino gradient applications.



---

**1 pcs**      **6.2744.350**      **PEEK pressure screw 14 mm, 3 pcs.**

With UNF 10/32 connector. For the connection of PEEK capillaries in the Absorber Module.

