





## Accessories for MATi 10

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery MATi 10

Qty.	Order no.	Description	
1 pcs	2.741.0010	<b>741 Magnetic Stirrer</b>	
		Magnetic stirrer for sample changers. The magnetic stirrer is located beneath the sample rack.	
2 pcs	2.800.0010	<b>800 Dosino</b>	
		Drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).	
1 pcs	2.814.0030	<b>814 USB Sample Processor (1T/0P)</b>	
		USB Sample Processor with one workstation for the automatic processing of routine samples in series with small to medium quantities. Up to two pumps (membrane or peristaltic) and three dosing devices for Liquid Handling tasks can be connected. Because of the multitude of application variants, rack, stirrer, titration head, Swing Head and sample vessels must be tailored to the application and ordered separately.	
1 pcs	2.901.0010	<b>901 Titrando</b>	
		The Karl Fischer Titrando for the modern titration laboratory In addition to the Karl Fischer mode (KFT), the high-end titrator is also equipped with the SET mode and can thus perform titrations to a predefined endpoint. "iTrode ready": As a result of the high-resolution, galvanic measuring interface – now also for the first time with additional digital measuring input for the intelligent "iTrodes" electrodes – it is additionally possible to perform pH measurements with a maximum of precision. Thanks to the four MSB connectors, up to four Dosinos can be connected to the 901 Titrando. These can be used not only for titration, but also for the automated addition of sample, standard and auxiliary solutions such as solubility promoters. Furthermore, you have all of the Liquid Handling commands at your disposal, which means that the aspiration of exhausted reagent and the addition of fresh solvent can also be automated.	

---

**24 pcs**    **6.1432.210**    **Sample beaker, clear glass, 75 mL, 1 piece**

Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



---

**1 pcs**    **6.1608.023**    **Amber glass bottle / 1000 mL / GL 45**

For Exchange Units. Bottle for auxiliary solutions.



---

**1 pcs**    **6.1805.030**    **FEP tubing / M6 / 150 cm**

With light and kink protection



---

**24 pcs**    **6.1903.020**    **Stirring bar / 16 mm**

Stirring Bar with magnetic core, PTFE coated, length 16 mm








---

**1 pcs**    **6.2041.340**    **Sample rack 24 x 75 mL**

Sample rack for 24 x 75 mL samples, one-row.



<p><b>1 pcs</b></p>	<p><b>6.2151.000</b></p> <p>Controller cable</p>	<p><b>Cable USB A – mini-DIN 8-pin</b></p>	
<p><b>1 pcs</b></p>	<p><b>6.3032.210</b></p> <p>Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p>	<p><b>Dosing Unit 10 mL</b></p>	
<p><b>1 pcs</b></p>	<p><b>6.3032.250</b></p> <p>Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p>	<p><b>Dosing Unit 50 mL</b></p>	
<p><b>1 pcs</b></p>	<p><b>6.5610.020</b></p> <p>Karl Fischer equipment for the 824 Easy Sample Changer.</p>	<p><b>Karl Fischer Equipment for 824</b></p>	
<p><b>1 pcs</b></p>	<p><b>6.6056.251</b></p> <p><b>Tiamo 2.5 Light CD: 1 license</b></p> <p>tiamo™ 2.5 Light computer program for controlling a titration system.</p> <ul style="list-style-type: none"> <li>• Up to two Metrohm instruments (Titrino, Titrand, etc.) can be connected; balances and other generic (i.e. non-Metrohm) instruments can be included without limitation</li> <li>• Graphical method editor with numerous templates</li> <li>• Layout manager for individual monitor interface</li> <li>• Professional database with reevaluation</li> <li>• High-performance report generator</li> <li>• Data export as PDF file, CSV, SLK</li> <li>• No parallel titration</li> <li>• No data export in XML format for LIMS</li> <li>• Dialog languages: German, English, French, Italian, Spanish, Czech, Portuguese, Polish, Russian, Slovakian, Japanese, Chinese, Traditional Chinese</li> </ul>	<p><b>Tiamo 2.5 Light CD: 1 license</b></p>	

## Optional accessories

Order no.	Description
<b>6.1459.400</b>	<b>Sample beaker, plastic (PE), 75 mL, 144 pieces</b>
	Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.

