






Accessories for 2.919.0130

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.919.0130

Qty.	Order no.	Description	
1 pcs	1.919.0020	919 IC Autosampler plus	
		The 919 IC Autosampler plus fulfills the requirements of laboratories with medium sample numbers. It enables automation of the full range of Metrohm IC instruments.	
1 pcs	6.1453.250	Sample beaker, plastic (PP), 250 mL, 1 piece	
		Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece	
1 pcs	6.1803.160	PTFE capillary 0.75 mm i.d. / 3 m	
		PEEK capillary for use with adsorber module.	
1 pcs	6.1815.010	Spiral band 0.5 m	
		For holding together different cables or tubing.	
2 pcs	6.1826.020	PVC pump tubing (blue/blue)	
		For online IC instruments and automation in voltammetry.	

1 pcs **6.1831.020** **PEEK capillary 0.75 mm i.d., 22 cm**
For all IC components.



1 pcs **6.2041.510** **Sample rack 56 x 11 mL+1 x 250 mL**
Sample rack for 56 x 11 mL samples with one rinsing beaker.



1 pcs **6.2151.000** **Cable USB A – mini-DIN 8-pin**
Controller cable



1 pcs **6.2621.080** **Capillary cutter**
For plastic capillaries. Used with IC instruments



1 pcs **6.2621.090** **1/2 in. wrench**
For 1/2 in. screws. For IC Sample Processors



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2621.130** **Hexagon key 2 mm**
2 mm.



1 pcs **6.2739.000** **Wrench**
For tightening connectors



1 pcs **6.2743.057** **Sample tubes 11 mL**
200 pieces. For IC Sample Processors and VA Autosampler



2 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



2 pcs **6.2744.160** **Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



1 pcs **6.2833.030** **Needle holder 1/8 in. for 863 Compact Autosampler**
Needle holder 1/8 in. for 863 Compact Autosampler



1 pcs **6.2846.000** **Sample needle made of zirconium oxide**

Sample needle for Sample Processors and Autosampler in ion chromatography and voltammetry.

