




Accessories for BA

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
RRDE.AUPT	Rotating ring disk electrode, Gold disk Platinum ring	
<p>The RRDE electrode tip consist of a 5 mm disc of gold with a concentric platinum ring at a distance of 375μm, leading to a theoretical collection efficiency of 24.9%.</p>		
RRDE.GCPT	Rotating ring disk electrode, Glassy Carbon disk Platinum ring	
<p>The RRDE electrode tip consist of a 5 mm disc of Glassy Carbon with a concentric platinum ring at a distance of 375μm, leading to a theoretical collection efficiency of 24.9%.</p>		
RRDE.PTPT	Rotating Ring Disk Electrode, Platinum disk Platinum ring	
<p>The RRDE electrode tip consist of a 5 mm disc of platinum with a concentric platinum ring at a distance of 375μm, leading to a theoretical collection efficiency of 24.9%.</p>		