



OMNIS Titrator Basic KFT

2.1001.6220

The OMNIS Titrator Basic KFT offers you the complete package for volumetric Karl Fischer titration. The package contains the OMNIS Titrator with magnetic stirrer for volumetric Karl Fischer titration, the Solvent Pump for manual reagent replacement, all necessary accessories for volumetric Karl Fischer titration, and one stand-alone OMNIS Software license. Titration can be carried out only with the internal buret of the OMNIS Titrator.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1001.6220

Order Number	Description	Quantity
1.1029.0010	Solvent Pump Solvent Pump is used with an Eco KF Titrator for reliably adding solvent and aspirating off titration vessel content without contact with the reagents.	1 PCS
6.01600.010	OMNIS Liquid Adapter Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	1 PCS
6.1546.030	Piston tongs For pistons for 807 Dosing Units or cylinder units OMNIS.	1 PCS

6.01601.000

Bottle cap, multi-use

1 PCS

Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.



6.2016.050

Support rod / 300 mm

1 PCS



6.2013.010

Clamping ring

1 PCS

For support rods with a diameter of 10 mm.



6.06000.000

OMNIS Software storage device

1 PCS

OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

License-based PC program for controlling complex OMNIS analysis systems.

- Sample-oriented operating procedures with numerous templates
- Graphical method editor with numerous templates
- Clearly organized and professional sample lists
- Management of instruments, work systems, sensors, solutions
- Dialog languages: German, English, Chinese

6.03004.210

Cylinder unit OMNIS special, 10 mL

1 PCS



Intelligent 10 mL cylinder unit for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Includes dosing tubing and antidiffusion tip.

6.02101.011

Measuring module analog

1 PCS



Analog measurement channel for one OMNIS Titrator or Titration Module for the connection of analog electrodes. Compatible with OMNIS software version 2.18 or higher.

6.06003.010

OMNIS Stand-Alone license

1 PCS

Enables stand-alone operation of the OMNIS software on a Windows™ computer.



Features:

- The license already includes one OMNIS instrument license.
- Must be activated via the Metrohm licensing portal.

6.2164.010

Power supply unit 100 - 240 V/24 V DC

1 PCS

Power supply unit for 915 KF Ti-Touch, 916 Ti-Touch, 917 Coulometer and all Eco Titrators / Eco Dosimat.



6.05609.020

KFT accessories for IONIC Titrator

1 PCS

Complete accessory set for volumetric Karl Fischer titration with the IONIC Titrator



Optional accessories

1.1042.0010

Dosino Module (cable 0.75 m)



The Dosino Module is a high-precision drive for the space-saving extension of an OMNIS titration system to include an additional intelligent buret for titrations, dosings, and other liquid handling tasks. With permanently attached cable (length 0.75 m).

1.1042.0020

Dosino Module (cable 2 m)



The Dosino Module is a high-precision drive for the space-saving extension of an OMNIS titration system to include an additional intelligent buret for titrations, dosings, and other liquid handling tasks. With permanently attached cable (length 2 m).
