



Bottle PP, 2 L, GL 45

6.01608.000

Bottle 2 L (PP) designed for OMNIS Reagent Organizer, with volume indication for eluent, regenerant, reagent, and rinsing solutions.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.01602.100

**Bottle cap GL 45, 2 x M6, UNF
10/32**

Bottle cap GL 45 for eluent, regenerant, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.



6.01834.010

**FEP tubing M6 with filter,
0.27 m**

FEP aspiration tubing with pre-installed particle filter (pore size 27 µm) for the 6.01608.000 bottle PP, 2 L, GL 45.



6.02837.030

CO2 Absorber (Tube)

CO2 Absorber (Tube) with a Luer adapter; protects reagents and eluents from particles and ensures consistent retention times for CO2-sensitive eluents.



6.02837.040

Filter tube

Filter tube with Luer adapter, protects eluents and reagents from particles.

