



OMNIS IC AS COOL Tray

6.02040.010

Cooled tray for accommodating sample racks in the OMNIS IC Autosampler. The tray maintains samples at a controlled low temperature, supporting the analysis of temperature-sensitive compounds. Rack positioning and handling correspond to the standard AS Tray, ensuring reliable and precise operation during automated sampling. The tray can be seamlessly integrated into the autosampler system and can be easily removed to allow convenient loading and preparation of sample racks directly on the laboratory bench before placing in the autosampler.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.02041.530

OMNIS sample rack 6 x 50 mL



Sample rack for holding 6 x 50 mL bottles in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.02041.540

OMNIS sample rack 12 x 22 mm



Sample rack for holding 12 x 22mm OD vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.02041.500

OMNIS sample rack 24 x 10 mL



Sample rack for holding 24 x 10 mL vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.02041.510

OMNIS sample rack 40 x 5 mL

Sample rack for holding 40 x 5 mL vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.02041.520

OMNIS sample rack 48 x 11 mm

Sample rack for holding 48 x 11 mm OD vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.02041.550

OMNIS sample rack, 2 x 200 mL

Sample rack for holding 2 x 200 mL bottles in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray