



OMNIS sample rack 48 x 11 mm

6.02041.520

Sample rack for holding 48 x 11 mm OD vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS/ASCOOL Tray.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.2743.100

**Vials PP, 0.7 mL, U-shape,
1000 pieces**



Maximize sample recovery with this 0.7 mL vial designed for low-volume analytical applications. The vial features a low-dead-volume insert design that enables efficient sample extraction while minimizing waste.

Quantity: 1000 pieces

6.2743.107

**Vials PP, 0.7 mL, U-shape, 100
pieces**



Maximize sample recovery with this 0.7 mL vial designed for low-volume analytical applications. The vial features a low-dead-volume insert design that enables efficient sample extraction while minimizing waste.

Quantity: 100 pieces

6.2743.110

**Vials PP, 0.3 mL, V-shape, 1000
pieces**



Maximize sample recovery with this 0.3 mL vial designed for low-volume analytical applications. The vial features a low-dead-volume insert design that enables efficient sample extraction while minimizing waste.

Quantity: 1000 pieces

6.2743.117

Vials PP, 0.3 ml, V-shape, 100 pieces




Maximize sample recovery with this 0.3 mL vial designed for low-volume analytical applications. The vial features a low-dead-volume insert design that enables efficient sample extraction while minimizing waste.

Quantity: 100 pieces

6.2743.120

Snap cap with red/white septum, 1000 pieces



Snap cap with septum for 11 mm vials. Provides reliable sealing and broad chemical compatibility for routine laboratory analysis

Quantity: 1000 pieces

6.2743.127

Snap cap with blue septum, 100 pieces



Snap cap with septum for 11 mm vials. Provides reliable sealing and broad chemical compatibility for routine laboratory analysis

Quantity: 100 pieces

6.7303.100

**Vials glass, 1.8 ml, U-shape,
100 pieces**



Maximize sample recovery with this universal 1.8 mL snap ring vial designed for low-volume analytical applications.

Quantity: 100 pieces

6.7303.110

**Snap cap with red septum,
100 pieces**



Snap cap with septum for 11 mm vials. Provides reliable sealing and broad chemical compatibility for routine laboratory analysis.

Quantity: 100 pieces