



Screw cap with septum, 12 pieces

6.02743.097

Screw cap with integrated septum designed for the secure storage and handling of samples, standards, and reagents. The septum enables automated sample access while minimizing contamination and evaporation, making it ideal for automated laboratory workflows. Its robust design ensures reliable sealing and dependable performance during sample preparation and analysis.

Quantity: 12 pieces

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.02743.030

Bottle, PP, 50 mL, 12 pieces



The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.

6.02743.040

Bottle PP, 200 mL



High-quality polystyrene (PS) sample vial designed for seamless use with the OMNIS IC Autosampler and for the safe handling and analysis of samples, reagents, and standards. Manufactured and packaged under controlled cleanroom conditions (ISO Class 7 according to DIN EN ISO 14644-1), the vial ensures low contamination levels and dependable performance in routine laboratory workflows.

Quantity: 1 piece
