



IC equipment: Inline Ultrafiltration 2 - pull mode

6.05330.010

Accessory set for Inline Ultrafiltration 2 - pull mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05330.010

Order Number	Description	Quantity
6.2744.070	Pressure screw short Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries	1 PCS
		
6.2621.030	Hexagon key 4 mm	1 PCS
		
6.1831.060	PEEK capillary 0.5 mm i.d., 1 m For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	1 PCS
		

6.1803.040



PTFE capillary 0.5 mm i.d. / 1 m

4 PCS

Capillary for sample handling in IC.

6.2744.030



Coupling olive/UNF 10/32, 4x

1 PCS

Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump

6.02057.030



Ultrafiltration cell 2 holder

1 PCS

Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.

6.02729.110



Ultrafiltration cell 2

1 PCS

For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.

6.1826.330

Pump tubing LFL (orange/white), 3 stoppers

1 PCS

For all IC instruments with peristaltic pump.



6.1826.390

Pump tubing LFL (yellow/yellow), 3 stoppers

1 PCS

For the sample solution in Inline Ultrafiltration.



6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces 1 PCS

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



6.2744.000

Pressure screw PVDF

1 PCS

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



6.2831.010

**Plastic tweezers**

1 PCS

For dialysis and filtration membranes

6.1803.100

**PTFE capillary, 0.5 mm ID / 40 cm**

1 PCS

Capillary for filtration and dialysis