



## Aspiration tip M8 with M10 screw nipple

6.1821.070

For aspirating reagent in combination with 803 Ti Stand and volumetric Karl Fischer cell. The aspiration tip fits tubing with an M8 thread. The attached tubing nipple has an M10 thread.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.1821.070

| Qt.   | Order no.  | Description                          |
|-------|------------|--------------------------------------|
| 1 PCS | 6.1808.230 | Adapter M6 / M8 inner                |
|       |            | Adapter with M6 and M8 inner threads |



1 PCS

6.1821.060

Aspiration tube with M6 thread

Spare part for 6.1821.070 aspiration tube

