



Al foil

6.2820.000

1,000 rondels with a diameter of 80 mm. Foil thickness: 0.010 mm. Together with 6.2037.040 foil holder for 6.1432.210 sample beaker and 6.2037.080 foil holder for 6.1459.300 sample beaker.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.2037.040

Foil holder



Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer

6.1432.210

Sample beaker, clear glass, 75 mL, 1 piece



Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece

6.2037.080

Foil holder



Foil holder for 120 mL PP beaker (6.1459.300)

6.1459.300

**Sample beaker, plastic (PP),
120 mL, 250 pieces**



Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces