




Rinsing/aspiration equipment for 869 Compact Sample changer or 862 Compact Titrosampler

6.5621.000

For usage of 869 Compact Sample Changer or 862 Compact Titrosampler in combination with a pump for rinsing/aspiration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5621.000

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1543.120

Aspiration tip / M8 thread

Aspiration tip if 6.1805.200 siphoning tubing is used



3 PCS

6.1805.020

FEP tubing / M6 / 52 cm

With light and kink protection



1 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.240

Distributor Rinsing / Aspiration for 869 Compact Sample Changer

Distributor with 3 x M6, 1 x M8 and 2 x 1/8 connections for the use of rinsing and / or aspiration equipment on the 869 Compact Sample Changer.



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.

