



QT-Sampler


BWT-840001283

The QT-Sampler is an accessory for transmission Raman measurements, enabling content uniformity testing of pharmaceutical tablets. Compatible with i-Raman NxG 785H.

Apertures and sample holders are sold separately.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery BWT-840001283

Order Number	Description	Quantity
BWT-840001177	Aperture Removal Tool Aperture removal tool for easy interchanging different apertures for optimal Raman spot size for transmission Raman with QT-Sampler. TOOL-APERTURE	1 PCS
		
BWT-840001153	QT-Sampler carrying case Lightweight black polycarbonate carrying case for QT-Sampler module. Watertight, crushproof, and dustproof with cube foam fitting. CCRM QT-Sampler	1 PCS
		

Optional accessories

BWT-840000940

QT-Sampler Optics and Sample Holder Set – 2 mm

Set of 2 mm excitation aperture, collection aperture, large round sample holder and large square sample holder for QT-Sampler.



BWT-840000941

QT-Sampler Optics and Sample Holder Set – 4 mm

Set of 4 mm excitation aperture, collection aperture, large round sample holder and large square sample holder for QT-Sampler.



BWT-840000942

QT-Sampler Optics and Sample Holder Set – 6 mm

Set of 6 mm excitation aperture, collection aperture, large round sample holder and large square sample holder for QT-Sampler.



BWT-840000934

Polystyrene Reference for QT-Sampler

Polystyrene reference for the QT-Sampler transmission Raman accessory.

