



# Flat cell


FLAT.CELL.S

The Autolab flat cell is a corrosion cell suitable for measurements on very large flat samples, up to 17 cm<sup>2</sup> in size. The cell is supplied with a reference electrode, a stainless steel counter electrode and a gas purging tube.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery FLAT.CELL.S


Qt.	Order no.	Description
1 PCS	FLAT.COUNTER	Flat cell steel counter electrode



1 PCS

FLAT.RING

Flat cell O-ring (set of 5)




---

1 PCS

FLAT.TOP

Flat cell PTFE top cover




---

1 PCS

FLAT.VESSEL

Flat cell vessel



---

1 PCS

6.0733.100

Ag/AgCl reference electrode

Silver / silver chloride reference electrode with  $c(\text{KCl}) = 3 \text{ mol/L}$  as reference electrolyte.

This reference electrode is suitable for aqueous applications and the SGJ 14/15 enables easy assembly.

