



703 Ti Stand with stirrer and pump

2.703.0010

納品範囲 2.703.0010

Qt.	Order no.	説明
1 PCS	1.703.0010	703 Ti Stand with stirrer and pump



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1446.060

Stopper / B-14/15 / M10



1 PCS

6.1543.110

Dosing tip / M8 thread



2 PCS

6.1602.105

Bottle attachment / GL 45 / green



1 PCS

6.1608.023

褐色ガラスボトル / 1000 mL / GL 45

交換ユニット用。補助溶液用ボトル。



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1609.010

Adsorber tube with tubing nipple



2 PCS

6.1801.000

PVC tubing / 4 mm / 7 mm / 1 m



1 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.050

Tubing adapter olive / M8 outer



1 PCS

6.1818.000

Aspiration tubing



1 PCS

6.1819.030

Tube



1 PCS

6.2001.010

Support stand with support rod



1 PCS

6.2013.010

Clamping ring



1 PCS

6.2023.020

Clip for SGJ 14/15



1 PCS

6.2108.100

Connection cable for 622/649/703/728, 0.5 m



1 PCS

6.2739.000

Wrench

