



705 UV Digester (220 V, 50 Hz)


2.705.0017

納品範囲 2.705.0017

Order Number	説明	数量
1.705.0017	705 UV Digester (220 V, 50 Hz)	1 PCS



1.705.0020	705 UV Digester wet end	1 PCS
------------	-------------------------	-------



6.2804.030	UV Mercury vapor lamp	1 PCS
------------	-----------------------	-------



6.1446.100

PTFE stopper for 12 mL UV quartz sample vessel 12 PCS

6.1446.110

**Thermometer stopper for UV quartz sample tube
12 mL** 1 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m 1 PCS

6.2041.040

Sample holder for 12 quartz digestion vessels, 12 mL 1 PCS



6.2401.010

Mercury thermometer

1 PCS



6.2414.000

Quartz sample vessel 12 mL


12 PCS



6.2745.000

UV protection plate

1 PCS



オプションアクセサリ

6.1446.150

**Thermometer stopper for
UV quartz sample tube 5
mL**



6.2041.050

**Sample holder for 16
quartz digestion vessels, 5
mL**



6.1446.140

**Stopper for 5 mL UV quartz
sample vessel**



6.2414.010

Quartz Sample Vessel 5 mL

