



## 786 Swing Head with transfer head right

2.786.0020

**納品範囲 2.786.0020**

| Order Number | 説明                                     | 数量    |
|--------------|----------------------------------------|-------|
| 1.786.0010   | 786 Swing Head with transfer head left | 1 PCS |



|            |                          |       |
|------------|--------------------------|-------|
| 6.1562.100 | Pipetting tubing / 10 mL | 1 PCS |
|------------|--------------------------|-------|



|            |                    |       |
|------------|--------------------|-------|
| 6.2621.120 | Hexagon key 1.5 mm | 1 PCS |
|------------|--------------------|-------|



6.1462.040

**Robotic arm with transfer head for 786 Swing  
Head, right swinging**

1 PCS



---

6.1823.010

**Guide shaft to pipetting tubing**

1 PCS

