



833 Advanced IC Liquid Handling Dialysis Unit

2.833.0040

納品範囲 2.833.0040

Order Number	説明	数量
6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m	7 PCS
		
6.1831.050	PEEK capillary 0.5 mm i.d., 40 cm	1 PCS
		
6.1831.060	PEEK capillary 0.5 mm i.d., 1 m	1 PCS
		

6.2057.010

Holder for dialysis cell

1 PCS



6.2128.180

Cable 762/830 - 752/753/754/833

1 PCS



6.2621.070

5 mm hex key for IC Sample Processors

1 PCS



6.2714.010

Dialysis membrane (cellulose acetate)

1 PCS



6.2729.100

Dialysis cell

1 PCS



6.2744.000

Pressure screw PVDF

1 PCS



6.2744.010

Pressure screw 5x

2 PCS



6.2744.034


カップリング チューブニップル/UNF 10/32、2個

2 PCS

圧力調整ネジとポンプチューブを接続。2個。ペリスタリックポンプ付きIC装置用。



6.2744.160

Nozzle/UNF 10/32 coupling with security device 4 PCS

6.2831.010

Plastic tweezers

1 PCS



オプションアクセサリー

6.1826.060
Pump tubing (orange
/yellow)



6.2714.030
Dialysis membrane
(polyamide)



6.2755.000
Tubing cartridge

