



843 Pump Station (membrane) - rinse/aspirate for Compact Sample Changer

2.843.0030

納品範囲 2.843.0030

Order Number	説明	数量
1.843.0020	843 Pump Station (membrane)	1 PCS
		
6.1543.120	Aspiration tip / M8 thread	1 PCS
		
6.1621.000	Container 10 L	2 PCS
		

6.1621.100

Cap for 10-L container

2 PCS



6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m

4 PCS



6.1805.200

PTFE tubing / M8 / 0.5 m

1 PCS



6.1805.420

FEP tubing / M6 / 48 cm

3 PCS



6.1808.240

Distributor Rinsing / Aspiration for 869 Compact 1 PCS
Sample Changer



6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m

2 PCS



6.1828.000

PVDF connection nipple

2 PCS



6.2740.020

Spray nozzle

3 PCS



オプションアクセサリ

6.2141.260

**Remote cable Compact
Sample Changer - Dosimat
plus and Titrino plus**



6.2141.300

**Remote cable 843 Pump
Station to Sample Processor**



6.2740.030

Rinsing nozzle

