



843 Pump Station (peristaltic)

2.843.0120

納品範囲 2.843.0120

Order Number	説明	数量
1.843.0120	843 Pump Station (peristaltic)	1 PCS



6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	4 PCS
------------	--------------------------------	----------



6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	2 PCS
------------	-------------------------------	----------



6.1820.050

M6 / Olive adapter

4
PCS


6.1826.100

Pump tubing 6.4 mm / 1 m

1
PCS


6.1828.000

PVDF connection nipple

2
PCS


オプションアクセサリ

6.2141.260

**Remote cable Compact
Sample Changer - Dosimat
plus and Titrino plus**



6.2141.300

**Remote cable 843 Pump
Station to Sample Processor**



6.2761.000

**Splash protection for 843
Pump Station (Peristaltic)**

