



843 Pump Station (peristaltic) - rinse/aspirate for Sample Processors

2.843.0150

納品範囲 2.843.0150

Order Number	説明	数量
1.843.0120	843 Pump Station (peristaltic)	1 PCS
		
6.1543.170	Aspiration tip / M8 thread	1 PCS
		
6.1621.000	Container 10 L	2 PCS
		

6.1621.100

Cap for 10-L container

2 PCS



6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m

4 PCS



6.1805.060

FEP tubing / M6 / 60 cm

3 PCS



6.1805.420

FEP tubing / M6 / 48 cm

3 PCS



6.1805.510

PTFE tubing / M8 / 60 cm

1 PCS



6.1808.170

Distributor for Sample Processors for rinsing and aspirating

1 PCS



6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m

2 PCS



6.1820.050

M6 / Olive adapter

4 PCS



6.1826.100

Pump tubing 6.4 mm / 1 m

1 PCS



6.1828.000

PVDF connection nipple

2 PCS



6.2141.300

Remote cable 843 Pump Station to Sample Processor

1 PCS



6.2740.020

Spray nozzle

3 PCS



オプションアクセサリー

6.2740.030

Rinsing nozzle

