



871 Advanced Bioscan

2.871.0010

納品範囲 2.871.0010

Order Number	説明	数量
1.871.0010	871 Advanced Bioscan	1 PCS



6.2156.010	Recorder cable for Bioscan	1 PCS
------------	----------------------------	----------



6.1816.060	Drain tubing for UV/VIS Compact IC and Bioscan	1 PCS
------------	--	----------



6.1825.210

PEEK sample loop 20 μ L

1
PCS


6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

1
PCS


6.2620.150

Pulsation absorber (metal-free)

1
PCS


6.2744.010

Pressure screw 5x

2
PCS


6.2744.014

Pressure screw 2x

1
PCS


6.2744.040

2 x UNF 10/32 coupling

1
PCS


オプションアクセサリ

6.1254.070

Silver flow cell for Bioscan



6.1254.130

Platinum Variocell for Bioscan



6.1254.140

Silver Variocell for Bioscan



6.1254.230

**Platinum working electrode
for Variocell**



6.2744.130

**Pressure screw for Bioscan
measuring cell**



6.1803.030

**PTFE capillary 0.5 mm i.d. /
3 m**

