

1 PCS

6.1414.060

マイクロ滴定容器上部 / ECO



1 PCS

6.1415.150

Titration vessel / 5-70 mL



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1448.010

Septum 12 mm, 5 pieces



1 PCS

6.1543.120

Aspiration tip / M8 thread



2 PCS

6.1805.010

FEP tubing / M6 / 13 cm



1 PCS

6.1805.200


PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.100

Distributor 3 x M6 / 1 x M8



1 PCS

6.1903.020

Stirring bar / 16 mm



2 PCS

6.2709.090

Stopper



1 PCS

6.2730.020

Septum stopper



5 PCS

6.2730.030

Stopper



2 PCS

6.2730.060

Screw nipple



2 PCS

6.2730.070

Screw nipple



1 PCS

6.2739.000

Wrench



2 PCS

6.2740.030

Rinsing nozzle



オプションアクセサリー

Order no.	説明
6.2730.080	Screw nipple

