



Transfer equipment, 0.5 - 2 mL, for 789/815

6.5619.000

納品範囲 6.5619.000

Order Number	説明	数量
6.1446.180	Stopper, PTFE, SGJ 29, for transfer set	1 PCS



6.1448.000	Septum stopper B-14/15	1 PCS
------------	------------------------	-------



6.1448.020	Septum 16 mm, 5 pieces	2 PCS
------------	------------------------	-------



6.1462.090

**Robotic arm with Luer-Lock adapter for 786
Swing Head, right swinging**

1 PCS



6.1464.320

カールフィッシャー滴定容器 / 80-250 mL / 電量滴定用 1 PCS



6.1562.110

Pipetting tubing / 2 mL

1 PCS



6.1608.090

Bottle / 250 mL / GL 45


1 PCS



6.2001.100

Base plate, transfer set

1 PCS



6.2013.010

クランプリング

1 PCS



直径10 mmの支持棒用。

6.2701.040

GL 18 screw cap with hole

1 PCS



6.2751.120

Safety shield Transfer set to 789/815

1 PCS



6.2816.100

Needle with vent

1 PCS



6.2055.110

**ボトル2本 / Dosino のためのボトルホルダー、もしくは
はDosinoホルダー** 1 PCS250 mL のボトル2本もしくは Dosino 2つのためのホル
ダー