



Autolab PGSTAT101

AUT101.S

納品範囲 AUT101.S

Qt.	Order no.	説明
1 PCS	CBL.USB	Standard USB cable



1 PCS

CABLE.PWR


Power cable



1 PCS

CELLCBL.101

Cell cable (CE, RE, S, WE, GND) for PGSTAT101



1 PCS

NOVA

Advanced software for electrochemical research

NOVA is the package designed to control all the Autolab instruments with USB interface. Designed by electrochemists for electrochemists and integrating over two decades of user experience and the latest .NET software technology, NOVA brings more power and more flexibility to your Autolab potentiostat /galvanostat. NOVA offers the following unique features: Powerful and flexible procedure editor; Clear overview of relevant real-time data; Powerful data analysis and plotting tools; Integrated control for external devices like Metrohm Liquid Handling devices; Download the latest version of NOVA




オプションアクセサリ

Order no.	説明
CABLE. MONITOR4	Monitor cable for PGSTAT101



AUT.DUMCELL. S	Autolab dummy cell
-------------------	--------------------




AUT101.CASE	AUT101.CASE
-------------	-------------



SDK	Software development kit
-----	--------------------------

Metro
Autolab



ALL.CLIP.BLACK Black Alligator Clamp



ALL.CLIP.RED Red Alligator Clamp

