





Autolab PGSTAT302N

AUT302N.S

This high end, high current potentiostat/galvanostat, with a compliance voltage of 30 V and a bandwidth of 1 MHz, combined with our FRA32M module, is specially designed for electrochemical impedance spectroscopy. The PGSTAT302N is the successor of the popular PGSTAT30. The maximum current is 2 A, the current range can be extended to 20 A with the BOOSTER20A, the current resolution is 30 fA at a current range of 10 nA.

納品範囲 AUT302N.S

Order Number	説明	数量
CELLCBL.30.RE.S	Cell cable	1 PCS
		
CELLCBL.30.WE.S	Cell cable	1 PCS
		
AUT.DUMCELL.S	Autolab dummy cell	1 PCS
		

CBL.USB

Standard USB cable

1 PCS



CABLE.PWR

Power cable

1 PCS



NOVA

Advanced software for electrochemical research

1 PCS



NOVA is the package designed to control all the Autolab instruments with USB interface. Designed by electrochemists for electrochemists and integrating over two decades of user experience and the latest .NET software technology, NOVA brings more power and more flexibility to your Autolab potentiostat /galvanostat. NOVA offers the following unique features: Powerful and flexible procedure editor; Clear overview of relevant real-time data; Powerful data analysis and plotting tools; Integrated control for external devices like Metrohm Liquid Handling devices; Download the latest version of NOVA

CABLE.MONITOR	Monitor cable for N series Autolab	1 PCS
---------------	------------------------------------	-------



CABLE.BNC.50	50 cm BNC cable	1 PCS
--------------	-----------------	-------



SDK	Software development kit	1 PCS
-----	--------------------------	-------

オプションアクセサリ

3500003600
20 A current booster



FRA32M.S
Electrochemical impedance
spectroscopy module



ADC10M.S
Dual-channel ultra-fast
sampling module



MBA.S
**Dual mode bipotentiostat
module**



ECN.S
**Electrochemical noise
module**



ECD.S
**Low current amplifier
module**



MUX.S
Multiplexer module



FI20.S
Analog filter and integrator module



pX1000.S
Voltage and pH measurement module



EQCM.S
Electrochemical Quartz
Crystal Microbalance
module




ECI10M.S
High frequency impedance
spectroscopy



BSTR10A.S
Booster 10A



IME663.S
Interface for Metrohm 663
VA Stand



ALL.CLIP.BLACK
Black Alligator Clamp



ALL.CLIP.RED
Red Alligator Clamp



SCAN250.S
True linear scan generator
module



BA.S
Dual mode bipotentiostat
module

