



Autolab PGSTAT100N

PGSTAT100N

納品範囲 PGSTAT100N

Qt.	Order no.	説明
1 PCS	CBL.USB	Standard USB cable



1 PCS

CABLE.PWR

Power cable



1 PCS

CABLE.
MONITOR


Monitor cable for N series Autolab



1 PCS

AUT.DUMCELL.
S

Autolab dummy cell



1 PCS


NOVA

Advanced software for electrochemical research



1 PCS

CELLCBL.30.WE. Cell cable
S



1 PCS

CELLCBL.30.RE.S Cell cable



1 PCS

CABLE.BNC.50 50 cm BNC cable



オプションアクセサリ

Order no.	説明
SDK	Software development kit

Metro Autolab

ADC10M.S Dual-channel ultra-fast sampling module



AUT302N.S Autolab PGSTAT302N



MBA.S Dual mode bipotentiostat module





FI20.S

Analog filter and integrator module



FRA32M.S

Electrochemical impedance spectroscopy module



ECD.S

Low current amplifier module



IME663.S

Interface for Metrohm 663 VA Stand



SCAN250.S

True linear scan generator module



ALL.CLIP.BLACK

Black Alligator Clamp



ALL.CLIP.RED

Red Alligator Clamp

