



OMNIS Sample Robot M – NIR

2.1074.0010

OMNIS Sample Robot M – NIR for the direct transition to fully automated spectroscopy. The system provides space for 3 sample racks for up to 231 vials with diameters of 8–28 mm.

The OMNIS Sample Robot M – NIR supports the analyses of both liquid and solid samples.

Phạm vi giao hàng 2.1074.0010

Qt.	Order no.	Description
1 PCS	1.1074.0010	OMNIS Sample Robot M – NIR

OMNIS Sample Robot M – NIR for the direct transition to fully automated spectroscopy. The system provides space for sample racks for up to 231 vials with diameters of 8–28 mm.



The OMNIS Sample Robot M – NIR supports the analyses of both liquid and solid samples.



Function license for utilization of the command parameters for the OMNIS Sample Robot M – NIR.



Phụ kiện tùy chọn

Order no.	Description	
6.02041.080	OMNIS sample rack 77 x 8 mm OMNIS sample rack for OMNIS Sample Robot – NIR suitable for 77 disposable vials with an optical path length of 8 mm for measurements in transmission (6.7402.240).	
6.02041.090	OMNIS sample rack 77 x 15 mm OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 77 disposable vials with a diameter of 15 mm for measurements in reflection (6.7402.110).	
6.02041.100	OMNIS sample rack 50 x 19 mm OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 50 disposable vials with a diameter of 19 mm for measurements in reflection (6.7402.120).	
6.02041.110	OMNIS sample rack 50 x 22 mm OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 50 disposable vials with a diameter of 22 mm for measurements in reflection (6.7402.130).	

6.02041.120

OMNIS sample rack 27 x 28 mm

OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 27 disposable vials with a diameter of 28 mm for measurements in reflection (6.7402.140).

