



Balance for 864 Robotic Balance Sample Processor, 115V

2.135.0205

The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 115V.

Phạm vi giao hàng 2.135.0205

Order Number	Description	Quantity
--------------	-------------	----------

1.135.0205

**Balance for 864 Robotic Balance Sample Processor,
115V**

1 PCS



The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 115V.

Phụ kiện tùy chọn

6.2067.100

Scale pan for sample rack
6.2068.010



The scale pan is designed for the 120 mL PP beaker (6.1459.300) and is placed into the sample rack of the 864 Balance Sample Processor.

6.2068.010

Sample rack, 20 x 120 mL,
(PP)



Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.

Plastic: Polypropylene (PP)

6.2711.080

Drip pan for 864 Robotic
Balance Sample Processor



The drip pan guards Robotic Balance Sample Processor and balance against contamination / corrosion by chemicals.