



Cylinder unit OMNIS 5 mL

6.03001.150

Intelligent 5 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.

Phạm vi giao hàng 6.03001.150

Order Number	Description	Quantity
6.01503.150	Cylinder unit OMNIS 5 mL without accessories Intelligent cylinder unit OMNIS 5 mL with four M6 ports for OMNIS Titrator, Titration Module, Dosing Module or OMNIS Dosing Drive Unit.	1 PCS
		
6.1446.030	Ball stopper / B-14/(15) For the movable holding of buret tips in titration vessels	1 PCS
		
6.1446.040	Threaded stopper / M6	2 PCS
		

6.1543.200


Titration tip / M6 thread

1 PCS

Titration tip with antidiffusion valve.

6.1805.100


FEP tubing / M6 / 40 cm

2 PCS

With light and kink protection

6.2008.030


Buret tip storage vessel for 807 Dosing Unit

1 PCS

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.

6.2709.070


Guiding sleeve

1 PCS

For fixation of buret and dosing tips in titration heads and electrode holders

6.2739.000

Wrench

1 PCS

For tightening connectors



Phụ kiện tùy chọn

6.05700.140

**Consumable Kit, cylinder unit
OMNIS**

FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.



6.05700.411

**Consumable Kit, cylinder unit
OMNIS 5 mL**

Cylinder unit OMNIS 5 mL complete with cylinder, piston, cylinder base, valve disk, centering tube, and cylinder housing, including programmed data chip.

