



Accessory kit: Bottle OMNIS IC, PP, 2 L

6.05000.010

OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.

Phạm vi giao hàng 6.05000.010

Order Number	Description	Quantity
6.01602.100	Bottle cap GL 45, 2 x M6, UNF 10/32 Bottle cap GL 45 for eluent, regenerant, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	1 PCS
6.01608.000	Bottle PP, 2 L, GL 45 Bottle 2 L (PP) designed for OMNIS Reagent Organizer, with volume indication for eluent, regenerant, reagent, and rinsing solutions.	1 PCS
6.01834.010	FEP tubing M6 with filter, 0.27 m FEP aspiration tubing with pre-installed particle filter (pore size 27 µm) for the 6.01608.000 bottle PP, 2 L, GL 45.	1 PCS

6.02050.000

Bottle holder, RO

1 PCS

Bottle holder for Eluent bottle PP 2 L on OMNIS IC Reagent Organizer.



6.02050.010

Bottle fixation RO

1 PCS

Bottle fixation for Eluent bottle PP 2 L on OMNIS IC Reagent Organizer.



Phụ kiện tùy chọn

6.05801.040

Volume sensor, ROLC kit



Upgrade kit to extend the OMNIS Reagent Organizer Liquid Control with an additional volume sensor.

Prerequisites:

– 2.1056.1200 OMNIS Reagent Organizer Liquid Control

Info: This upgrade must be installed by an authorized Metrohm service technician.

2.1056.1200

OMNIS Reagent Organizer Liquid Control



The OMNIS Reagent Organizer Liquid Control serves as Reagent Organizer and provides the illumination of the interior in the OMNIS IC System. In addition, it offers a software-integrated leak detection.

The OMNIS Reagent Organizer Liquid Control is always located on the topmost OMNIS Extension Module and also forms the lid of the OMNIS IC System. A leak sensor is installed at the drain of the OMNIS Reagent Organizer Liquid Control. The lighting integrated into the OMNIS Reagent Organizer Liquid Control illuminates the interior of the OMNIS IC System as soon as at least one door is opened.

The OMNIS Reagent Organizer Liquid Control offers space for up to 6 bottles and can be equipped with additional accessories, as for example bottle holder, bottle fixation, volume sensor, magnetic stirrer or waste container with external canister sensor.

6.02061.000

Reagent Organizer



Organizer for reagent and eluent bottles for OMNIS IC instruments.