



Accessory kit: Vario/Flex ONE

6.5000.010

Set of accessories for the ONE Module (IC pump with injector) for Professional IC Vario and Compact IC Flex instruments.

Phạm vi giao hàng 6.5000.010

Order Number	Description	Quantity
6.02744.050	Adapter Luer inner / M6 outer Adapter for syringe for pre-rinsing the aspiration filter (6.2821.090). The adapter can be screwed onto the aspiration filter. A syringe can be used to rinse the aspiration filter with ultrapure water or eluent.	1 PCS
6.1602.160	Bottle attachment GL 45 for eluents and reagents For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	1 PCS
6.1608.070	Eluent bottle / 2 L / GL 45	1 PCS



6.1609.000

Adsorbing tube, large and bent

1 PCS

For filling with adsorber material.



6.1803.040

PTFE capillary 0.5 mm i.d. / 1 m

1 PCS

Capillary for sample handling in IC.



6.2023.020

Clip for SGJ 14/15

1 PCS

Clip for SGJ 14/16



6.2744.014

Pressure screw 2x

1 PCS

With UNF 10/32 connection. For the connection of PEEK capillaries



6.2744.040

2 x UNF 10/32 coupling

1 PCS

For connecting 1/16 in. capillaries. For IC instruments



6.2744.070

Pressure screw short

1 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.2744.090

Pressure screw long

1 PCS

Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



6.2744.210

Tubing adapter for aspiration filter (ProfIC)

1 PCS

For Professional IC instruments



6.2816.020

Syringe 10 mL with Luer connection

1 PCS

For various applications in IC and VA



6.2816.040

Purge tubing

1 PCS

With PTFE tubing and Luer connection. For syringes. For aspirating eluents.



6.2821.090

Aspiration filter

1 PCS

Pore size: 27 µm. Set of 5 pieces. For aspiration tubing (6.1834.000, 6.01834.010) and filter tubes (6.1821.040, 6.1821.050).



6.2821.130

Spare filter for inline filter

1 PCS

Spare filter for inline filter.

