




Rinsing/aspiration equipment for 869 Compact Sample changer or 862 Compact Titrosampler

6.5621.000

For usage of 869 Compact Sample Changer or 862 Compact Titrosampler in combination with a pump for rinsing/aspiration.

Phạm vi giao hàng 6.5621.000

Qt.	Order no.	Description
1 PCS	6.1543.120	Aspiration tip / M8 thread
		Aspiration tip if 6.1805.200 siphoning tubing is used



3 PCS 6.1805.020 FEP tubing / M6 / 52 cm
With light and kink protection



1 PCS 6.1805.200 PTFE tubing / M8 / 0.5 m



1 PCS 6.1808.240 Distributor Rinsing / Aspiration for 869 Compact
Sample Changer

Distributor with 3 x M6, 1 x M8 and 2 x 1/8 connections for the use of rinsing and / or aspiration equipment on the 869 Compact Sample Changer.



3 PCS 6.2740.020 Spray nozzle

For the fine-spraying of the rinsing solution.

