



Consumable Kit IC

6.5703.010

Assembly of important expendable items for ion chromatography.

Phạm vi giao hàng 6.5703.010

Order Number	Description	Quantity
6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m For all IC instruments	1 PCS
		
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	1 PCS
		
6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m For all IC components.	1 PCS
		

6.1831.030


PEEK capillary 0.75 mm i.d. / 3 m

1 PCS

For all IC components.

6.1831.180


PEEK Capillary 0.5 mm i.d., 3 m

1 PCS

6.2251.000


Colored sleeves for capillaries

1 PCS

Colored pieces of heat shrink tubing for capillary coding.
Three pieces each of five different colors.

6.2621.080


Capillary cutter

1 PCS

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.

6.2744.010

Pressure screw 5x

1 PCS

With UNF 10/32 connection. For the connection of PEEK capillaries



6.2821.090

Aspiration filter

1 PCS

Pore size: 27 µm. Set of 5 pieces. For aspiration tubing (6.1834.000, 6.01834.010) and filter tubes (6.1821.040, 6.1821.050).




6.2821.130

Spare filter for inline filter

1 PCS

Spare filter for inline filter.



6.2744.090

Pressure screw long

1 PCS

Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



6.2744.070

Pressure screw short

1 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries

