



Vials glass, 1.8 ml, U-shape, 100 pieces

6.7303.100

Maximize sample recovery with this universal 1.8 mL snap ring vial designed for low-volume analytical applications.

Quantity: 100 pieces

Phụ kiện tùy chọn

6.02041.520

OMNIS sample rack 48 x 11 mm



Sample rack for holding 48 x 11 mm OD vials in the OMNIS IC Autosampler. Rack geometry is matched to vial dimensions to ensure precise positioning and stable handling during automated sampling. The racks can be seamlessly mounted on the AS /ASCOOL Tray.

6.2041.200

Sample rack, 48 x 11 mL, (Al)




Sample rack for 889 IC Sample Center, suitable for 48 sample vessels. The following sample vessels can be used: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117.

Material: Aluminum (Al)

6.7303.110

Snap cap with red septum, 100 pieces




Snap cap with septum for 11 mm vials. Provides reliable sealing and broad chemical compatibility for routine laboratory analysis.

Quantity: 100 pieces

6.7303.310

**Sample rack 112 positions for
liquid samples**

For liquid samples.



6.7303.000

**Kit for liquid samples for
MMS 5000, 112 positions**

Kit for liquid samples for MMS 5000 with rack for sample vessels,
100 μ L syringe with replaceable needle and 100 snap-ring bottles.

