



826 pH mobile, 帶便攜箱

2.826.0110

便携式 pH 计，带 IrDA 接口。带便携箱和用于校正的多种附件，包括 Primatrode 6.0228.020。

Scope of delivery 2.826.0110

Order Number	Description	Quantity
1.826.0010	826 pH mobile 便携式 pH 计，带 IrDA 接口，包括背带和电池。	1 PCS
6.2050.000	pH计的背带 用于704，826型pH计的背带。	1 PCS
6.2716.040	826 pH计的机盒 826 pH计及配件的机盒	1 PCS



6.1236.050

套管, SGJ 14/12 毫米

1 PCS

套管, SGJ 14/12 毫米, 聚乙烯



6.1446.000

SGJ 塞子 / B-14/ (15)

1 PCS



6.1613.020

试剂瓶 / 25毫升 / 带pH 7标签 (空瓶)

1 PCS



6.1613.030

试剂瓶 / 25毫升 / 带pH 4标签 (空瓶)

1 PCS



6.1614.000

洗瓶 / 250毫升

1 PCS



6.1614.030

塑料瓶 / 50毫升

3 PCS



塑料瓶

6.2133.000

5号碱性电池，1.5V，4只

1 PCS



每套包含4节1.5伏的5号碱性电池。

6.2307.230

缓冲溶液，pH 4，7与9

1 PCS



混合缓冲液， pH 4.00/7.00/9.00 (25°C)， 装于单次使用的小袋中，无色，盒装，每盒3 x 10 x 30 毫升

6.2308.050

电解液, 3 mol/L KCl (50 mL)

1 PCS

电解液c (KCl) = 3 mol/L, 50 mL (用于Ag/AgCl参比系统)



6.2621.140

六角扳手 (Hexagon key), 2.5毫米

1 PCS



6.2717.000

PP 杯, 100 mL

1 PCS

PP 聚丙烯制成的量杯, 100 mL。



Optional accessories

6.2325.000
pHit kit

pH电极的维护套件该套件包括:- 50 mL 清洗液- 50 mL 3M KCl
溶液- 50 mL 储藏液- 2 储存容器- 使用说明

