



827 pH lab; 240 V, AUS 带 Primatrode

2.827.0117

实验室 pH 计，带 IrDa 接口和 Primatrode 6.0228.010

Scope of delivery 2.827.0117

Order Number	Description	Quantity
6.2161.040	插入式电源单元, AUS 240V/6V DC 电源, 用于827型 pH lab	1 PCS



6.2013.010	定位环	1 PCS
------------	-----	-------



用于直径为 10 mm 的支架杆。

6.2016.050	支撑杆 / 300毫米	1 PCS
------------	-------------	-------



6.2021.020

电极架

1 PCS

电极架，可用于4只电极与2只滴定头



6.2621.070

内六角扳手，5 mm，用于 IC Sample Processor

1 PCS



Optional accessories

6.2307.230
缓冲溶液, pH 4, 7与9

混合缓冲液, pH 4.00/7.00/9.00 (25°C), 装于单次使用的小袋中, 无色, 盒装, 每盒3 x 10 x 30 毫升



6.2325.000
pHit kit

pH电极的维护套件该套件包括:- 50 mL 清洗液- 50 mL 3M KCl 溶液- 50 mL 储藏液- 2 储存容器- 使用说明

