

807 Dosing Unit 2 mL

6.3032.120

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.3032.120

Order Number	Description	Quantity
6.1446.030	**Ball stopper / B-14/(15)** For the movable holding of buret tips in titration vessels	1 PCS
6.1446.040	**Threaded stopper / M6**	1 PCS
6.1543.200	**Titration tip / M6 thread** Titration tip with antidiffusion valve.	1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

1 PCS

With light and kink protection



6.2008.030

Buret tip storage vessel for 807 Dosing Unit

1 PCS

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



6.2803.000

Grease (1 oz.)

0 PCS

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



6.2739.000

Wrench

1 PCS

For tightening connectors



6.1543.060

Tip / M6 thread

1 PCS



Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.

6.1619.000

Adsorber tube for 807 Dosing Unit

1 PCS



6.1829.010

FEP aspiration tube M6, 0.25 m

1 PCS



With M 6 thread. For Dosing Unit

6.2008.050

Storage vessel holder for 807 Dosing Unit

1 PCS



Together with the 6.2008.030 buret tip storage vessel, it provides a holder for the buret tip on the 807 Dosing Unit.

6.2244.020

Cards for intelligent Exchange Unit 1 PCS



Set of 10 pieces, different colors.

6.2803.010

Grease (2 g) 1 PCS



Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

Optional accessories

6.1608.030
Round glass bottle / 1000 mL / GL 45



6.1608.040
Bottle PE / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



6.1608.050
Drying bottle / 100 mL / GL 45



6.1608.060
Bottle / 100 mL




6.1618.000
Thread adapter 32 mm to GL 45

Adapter for bottles with 32 mm thread (Fluka, Riedel de-Haen 500 mL...) to GL45



6.1618.010
Thread adapter 28 mm to GL 45

Adapter for bottles with 28 mm thread (Fisher...) to GL45



6.1618.020
Thread adapter S 40 to GL 45

Adapter for bottles with S40 thread (Merck...) to GL45.



6.1618.050

Thread adapter 40 mm to GL 45



Adapter for bottles with 40 mm thread to GL45.

6.1618.060

Thread adapter GL 45 to GL 45



Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange

---

6.1608.023

Amber glass bottle / 1000 mL / GL 45



For Exchange Units. Bottle for auxiliary solutions.

6.1808.280
Adapter 807 Dosing Unit Port 4, M6 inner

Adapter for using Port 4 of the 807 dosing Unit for connection with an FEP tubing / M6.



6.1574.120
Cylinder element glass 2 mL

Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 2 mL.